Package 'readNSx'

October 24, 2025

Title Read 'Blackrock-Microsystems' Files ('NEV', 'NSx')

Version 0.0.6

Description Loads 'Blackrock' https://blackrockneurotech.com neural signal data files into the memory, provides utility tools to extract the data into common formats such as plain-text 'tsv' and 'HDF5'.

License MPL-2.0 | file LICENSE

Language en-US

Encoding UTF-8

RoxygenNote 7.3.2

Imports data.table, fastmap, hdf5r, jsonlite, R6

LinkingTo cpp11

Suggests bit64, tools, testthat (>= 3.0.0), knitr, rmarkdown, spelling

SystemRequirements HDF5 (>= 1.8.13), little-endian platform

NeedsCompilation yes

Copyright For the readNSx package: Zhengjia Wang.

URL http://dipterix.org/readNSx/

BugReports https://github.com/dipterix/readNSx/issues

Config/testthat/edition 3

VignetteBuilder knitr

Author Zhengjia Wang [aut, cre]

Maintainer Zhengjia Wang <dipterix.wang@gmail.com>

Repository CRAN

Date/Publication 2025-10-24 05:20:02 UTC

2 get_channel

Contents

Index	10
	read_bci2000
	import_nsp
	get_specification
	get_nsx
	get_nsp
	get_nev
	get_file_type
	get_event
	get_channel

get_channel Get channel data

Description

Obtain channel information and data from given prefix and channel ID.

Usage

```
get_channel(x, channel_id)
```

Arguments

x path prefix specified in import_nsp, or 'nev/nsx' object

channel_id integer channel number. Please be aware that channel number, channel ID, electrode ID refer to the same concept in 'Blackrock' 'NEV' specifications. Electrodes are not physical metals, they refer to channels for historical reasons.

Value

A list containing channel data and meta information, along with the enclosing 'NSx' information; for invalid channel ID, this function returns NULL

get_event 3

get_event

Get event data packets from 'NEV'

Description

Get event data packets from 'NEV'

Usage

```
get_event(x, event_type, ...)
```

Arguments

Χ path prefix (see import_nsp), or 'nev/nsx' object event type to load, common event types are event_type 'digital_inputs' packet identifier 0 'spike' packet identifier 1 to 10000 as of version 3.0 'recording' packet identifier 65529 as of version 3.0, available after version 'configuration' packet identifier 65530 as of version 3.0, available after version 3.0 'log' packet identifier 65531 as of version 3.0, available after version 3.0 'button_trigger' packet identifier 65532 as of version 3.0, available after version 3.0 'tracking' packet identifier 65533 as of version 3.0, available after version 'video_sync' packet identifier 65534 as of version 3.0, available after version 'comment' packet identifier 65535 as of version 3.0, available after version 2.3 pass to other methods

Value

A data frame of corresponding event type, or NULL if event is not found or invalid

get_nev

get_file_type

Get 'Blackrock' file type

Description

Reads the first 10 bytes containing file type and version information.

Usage

```
get_file_type(path)
```

Arguments

path

path to the 'Blackrock' '.nev' or '.nsx' file, or a binary connection.

Value

A list containing file information, including file type, version information, and normalized absolute path.

get_nev

Load 'NEV' information from path prefix

Description

Load 'NEV' information from path prefix

Usage

```
get_nev(x, ...)
```

Arguments

```
x path prefix specified in import_nsp, or 'nev/nsx' object
... reserved for future use
```

Value

'NEV' header information if x is valid, otherwise NULL. See Section "'NEV' Data" in get_specification

get_nsp 5

get_nsp

Get a collection list containing 'NEV' and 'NSx' headers

Description

Get a collection list containing 'NEV' and 'NSx' headers

Usage

```
get_nsp(x)
```

Arguments

Х

path prefix specified in import_nsp, or 'nev/nsx' object

Value

A list containing 'nev' and imported 'nsx' headers, see import_nsp for details

get_nsx

Load 'NSx' information from path prefix

Description

Load 'NSx' information from path prefix

Usage

```
get_nsx(x, which, ...)
```

Arguments

```
x path prefix specified in import_nsp, or 'nev/nsx' object
which which 'NSx' to load, for example, which=3 loads 'ns3' headers
... reserved for future use
```

Value

'NSx' header information if data is found, otherwise returns NULL. See Section "'NSx' Data" in $get_specification$

6 get_specification

<pre>get_specification Get '.nev' or</pre>	'nsx	' specification
--	------	-----------------

Description

```
Get '.nev' or 'nsx' specification
```

Usage

```
get_specification(version, type = c("nev", "nsx"))
```

Arguments

version either character string or a vector of two integers; for example, "2.2", "2.3",

c(3, 0). Currently only these three versions are supported since I was unable to find any other versions. Please file an issue ticket if others versions are desired.

type file type; choices are 'nev' or 'nsx'

Value

The file specification as a list. The specification usually contains three sections: basic header (fixed length), extended header (dictionary-style), and data packets (data stream). The specification is used to parse the data files.

'NEV' Data

A 'NEV' object consists of three sections:

Section 1 contains basic information such as the time-origin of all the time-stamps, the time-stamp sampling frequency, data packets sizes.

Section 2 is extended header containing the configurations of channels, digital signals, etc. For any data packets in section 3, there should be at least one table in this section describing the settings.

section 3 is a collection of event packets such as digital signal inputs (most likely to be used at version 2.2 or by 'Ripple'), spike waveform, comments (sometimes storing epoch information), etc.

Please be aware that while most common entries can be found across different file versions, some entries are version-specific. If you are making your script general, you need to be very careful handling these differences. For more information, please search for the data specification manual from the 'Blackrock-Microsystems' website.

'NSx' Data

A 'NSx' file refers to the data files ending with 'ns1' through 'ns9'. Common types are 'ns2' (sampling at 1000 Hz), 'ns3' (sampling at 2000 Hz), and 'ns5' (sampling at 30,000 Hz).

A 'NSx' file also consists of three sections. Section 1 contains basic information such as the timeorigin of all the time-stamps, sampling frequencies, and channel counts within the file. Please be careful that item time_resolution_timestamp is not the sampling frequency for signals. This item import_nsp 7

is the sampling frequency for time-stamp. To obtain the signal sample rate, divided time_resolution_timestamp by period. For example, 'ns3' usually has time-stamp resolution 30,000 and period=15, hence the signal sample rate is 30000/15=2000Hz.

Section 2 usually contains one and only one channel table of which the number of rows should coincide with number of channels from section 1. Other information such as channel labels, physical connectors, pins, units, filter settings, digital-to-analog conversion are also included. Since readNSx always attempts to convert signals in 'volts' or 'milli-volts' to 'micro-volts', the 'units' column might be different to what's actual recorded in the 'NSx' file headers.

Section 3 contains partitions of continuous recording. When imported/loaded from readNSx, the digital signals are always converted to analog signals with 'micro-volts' unit. Please use get_channel to get the channel data.

Examples

```
get_specification(version = c(2,3), type = "nev")
get_specification(version = "3.0", type = "nsx")
```

import_nsp

Import signal data from 'Blackrock-Microsystems' data files

Description

Please use import_nsp to import 'NEV' and 'NSx' files.

Usage

```
import_nsp(
  path,
  prefix = NULL,
  exclude_events = "spike",
  exclude_nsx = NULL,
  verbose = TRUE,
  partition_prefix = "/part"
)
```

Arguments

```
path path to 'NEV' or 'NSx' files

prefix path prefix to save data files into

exclude_events exclude one or more 'NEV' data events, choices are 'spike', 'video_sync',
    'digital_inputs', 'tracking', 'button_trigger', 'comment', 'configuration';
    default is 'spike' since the spike data takes very long time to parse, and most
    users might still use their own offline spike-sorting algorithms.
```

8 import_nsp

```
exclude_nsx excluded 'NSx' types, integer vectors from 1 to 9; for example, exclude_nsx=c(1,2) will prevent 'ns1' and 'ns2' from being imported

verbose logical or a progress object: when logical, verbose indicates whether to print the progress as messages; when verbose is a progress object, this object must have inc method; examples are Progress in the shiny package, progress2 from dipsaus, or shiny_progress from shidashi

partition_prefix additional prefix to the data partition; default is "/part"
```

Value

A list of configurations, see get_specification for what's contained.

Examples

```
# Please get your own sample data first. This package does not
# provide sample data for privacy and license concerns :)
if(interactive() && file.exists("sampledata.nev")) {
library(readNSx)
# ---- Import for the first time ------
import_nsp(
 path = "sampledata.nev",
 prefix = file.path(
    "~/BIDSRoot/MyDataSet/sub-YAB/ses-008/ieeg/",
    "sub-YAB_ses-008_task-congruency_acq-NSP1_run-01"
 ),
 exclude_events = "spike", partition_prefix = "/part"
# ---- Load header information ------
prefix <- "sub-YAB_ses-008_task-congruency_acq-NSP1_run-01"</pre>
nev <- get_nev(prefix)</pre>
ns3 <- get_nsx(prefix, which = 3)</pre>
# get nev from nsx, or nsx from nev
get_nev(ns3)
get_nsx(nev, which = 5)
# ---- Load channel data
result <- get_channel(prefix, channel_id = 10)</pre>
channel_signal <- result$channel_detail$part1$data</pre>
channel_signal[]
}
```

read_bci2000

read_bci2000

Read 'BCI2000' recording data

Description

```
Read 'BCI2000' recording data
```

Usage

```
read_bci2000_header(file)
read_bci2000(file)
```

Arguments

file

path to the recording data

Value

Parsed signal data

Examples

```
# Package comes with sample data
file <- system.file("samples", "bci2000_sample.dat", package = "readNSx")
result <- read_bci2000(file)

print(result)

# Notive: v1.0 and v1.1 are different, but all in `Source` section
# sample rate
result$parameters$Source$SamplingRate$value

# Signal data 64 channels x 500 time-points
dim(result$signals)</pre>
```

Index

```
get_channel, 2, 7
get_event, 3
get_file_type, 4
get_nev, 4
get_nsp, 5
get_nsx, 5
get_specification, 4, 5, 6, 8
import_nsp, 2-5, 7
read_bci2000, 9
read_bci2000_header (read_bci2000), 9
```